



Connector Ring

Product Images



The Connector Ring is a sterile device to support active filtration, cell detachment, cell lysis or cell stimulation in combination with a [pluriStrainer](#) and a 50ml standard (conical) centrifuge tube. pluriStrainer, Connector Ring and a 50ml tube can be combined to a unit. This unit allows for a controlled flow through the sieve. When the Luer-lock of the Connector Ring is closed, sample or liquid remains on the sieve. An open Luer-lock starts the separation or the filtration process.

Additional Information

| | |
|-------|--|
| 大小 | 25 pcs. |
| 交付时间 | 1-2 |
| 变化 | 25 pcs, sterile |
| 储存条件 | 室温 |
| 规范性声明 | For research use only. Not for use in diagnostic procedures, |

