



## S-pluriBead® Maxi Reagent Kit

### Product Images



Description

Maxi S-pluriBead Reagent Kit contains 100ml 10xWash buffer, stabilization buffer, incubation buffer,detachment activation buffer, 10x 1,8ml detachment buffer concentrate, 10 S-pluriStrainer, 10 Funnel, 10 connector rings

Additional Information

大小	10
交付时间	1-2
品种	人类
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	Room Temperature; Buffer 4-8 °C; Buffer D -20°C
规范性声明	For research use only. Not for use in diagnostic procedures,

