



## CD11b M-pluriBead® anti-ms

### Product Images



Description

CD11b pluriBead® suspension for the positive cell isolation of CD11b+ target cells from any murine sample material. (2ml)

Additional Information

|       |                                                                               |
|-------|-------------------------------------------------------------------------------|
| 大小    | 2 ml                                                                          |
| 品种    | 小鼠(鼠类)                                                                        |
| 细胞类型  | 单核细胞, 粒细胞, 巨噬细胞, NK细胞                                                         |
| 样品料   | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| 纯度    | >90%                                                                          |
| 储存条件  | 4-8 °C                                                                        |
| 规范性声明 | For research use only. Not for use in diagnostic procedures,                  |

