



CD326 M-pluriBead® anti-ms

Product Images



Description

CD326 pluriBead® suspension for the positive cell isolation of CD326+ target cells from any murine sample material. (2ml)

Additional Information

大小	2 ml
品种	小鼠(鼠类)
细胞类型	循环肿瘤细胞, CTC (循环肿瘤细胞), 上皮细胞
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures,

