



CD11b S-pluriBead® anti-hu

Product Images



Description

CD11b pluriBead® suspension for the positive cell isolation of CD11b+ target cells from any human sample material. (2ml)

Additional Information

大小	2 ml
品种	人类
细胞类型	单核细胞, 粒细胞, 巨噬细胞, NK细胞
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures,

