




## 10x PBS Solution

### Product Images



Description

PBS (phosphate buffered saline) is a balanced salt solution used for a variety of cell culture applications, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting and preparing reagents. PBS is formulated without calcium and magnesium for rinsing chelators from the culture before cell dissociation.

This PBS is manufactured **without:**

- Calcium
- Magnesium
- Phenol Red
- Protease free

Components

- NaCl (1368,9 mMol/L)
- Na 2 HPO 4 •2H 2 O (100mMol/L)
- KH 2 PO 4 (19,8 mMol/L)
- KCl (26,8 mMol/L)

Additional Information

交付时间	1-2
感兴趣的领域	细胞生物学
不育症	无菌型
相位范围	pH 6.8
运输条件	室温
储存条件	室温
规范性声明	For research use only. Not for use in diagnostic procedures,

Product Options

尺寸:	100 ml
	500 ml

