



ROLLER 10 basic

Product Images



描述

Roller shaker with 10 rolls provides smooth rocking and rolling action at a fixed speed. Used for mixing blood samples, solid and liquid suspensions and viscous samples. The unit can be used with all conventional tubes and cylindrical bottles. Easily removable attachments provide for quick cleaning in case of sample spills.

- Usable at 4 to 50 °C in incubators
- Suitable for continuous operation
- Removable rolls enable shaking of larger tubes
- Easy operation with touch keypad
- Shaking diameter: 24,5 mm (height)
- Speed range: 30 rpm (fixed)

Additional Information

| | |
|------|-----|
| 交付时间 | 1-2 |
| 运输条件 | 室温 |
| 储存条件 | 室温 |

