



## ÜberStrainer 60 µm

### Product Images

---



Description

The ÜberStrainer is more than a cell strainer - ÜberStrainer is a sterile Sample Preparation Device to separate, isolate, concentrate and process particle out of complex and/ or diluted suspensions. The air tight screw cap with a Luer Lock connector makes it possible to filtrate an unlimited amount of liquid. The strainer can be easily separated from the rest of the device, fits into a 50 ml tube and into the 20 mm tissue plate cavity.

Additional Information

网眼尺寸	60 µm
大小	10 pcs
布料	PET(聚乙二醇酯)
外壳材料	PP（聚丙烯）
Mesh Type	woven
不育症	无菌型
储存条件	室温
规范性声明	For research use only. Not for use in diagnostic procedures.,

