



Syringe Strainer 70 µm

Product Images



Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

Additional Information

交付时间	1-2
网眼尺寸	70 μm
大小	12 pcs.
流量控制	Yes
外壳材料	LD-PE (低密度聚乙烯)
Mesh Type	woven
不育症	无菌型
运输条件	室温
储存条件	室温
规范性声明	For research use only. Not for use in diagnostic procedures.

Product Options

Variation:	12 pcs sterile
	12 pcs, non-sterile

