



## Connector Ring

### Product Images

---



## Description

---

The Connector Ring is a sterile device to support active filtration, cell detachment, cell lysis or cell stimulation in combination with a [pluriStrainer](#) and a 50ml standard (conical) centrifuge tube. pluriStrainer, Connector Ring and a 50ml tube can be combined to a unit. This unit allows for a controlled flow through the sieve. When the Luer-lock of the Connector Ring is closed, sample or liquid remains on the sieve. An open Luer-lock starts the separation or the filtration process.

## Additional Information

---

大小	25 pcs.
交付时间	1-2
变化	25 pcs, sterile
储存条件	Room Temperature
规范性声明	For research use only. Not for use in diagnostic procedures,

