



universal S-pluriBead® anti-rabbit

Product Images



Description

pluriBead universal are developed for the use with external antibodies. They are coated with a polyclonal antibody produced from the serum of a host animal. (2ml)

Additional Information

大小	2 ml
品种	Rabbit
细胞类型	Any
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures,

