



## CD43 S-pluriBead® anti-hu

### Product Images



Description

CD43 pluriBead® suspension for the positive cell isolation of CD43+ target cells from any human sample material. (2ml)

Additional Information

|       |                                                                               |
|-------|-------------------------------------------------------------------------------|
| 大小    | 2 ml                                                                          |
| 品种    | Human                                                                         |
| 细胞类型  | Monocytes, T cells, Granulocytes                                              |
| 样品料   | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| 纯度    | >90%                                                                          |
| 储存条件  | 4-8 °C                                                                        |
| 规范性声明 | For research use only. Not for use in diagnostic procedures,                  |

