



## CD4 S-pluriBead® anti-hu

### Product Images



Description

CD4 pluriBead® suspension for the positive cell isolation of CD4+ target cells from any human sample material. (2ml)

Additional Information

|       |   |
|-------|---|
| 大小    | 2 ml  |
| 品种    | Human   |
| 细胞类型  | Monocytes (subset), Helper T cells  |
| 样品料   | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| 纯度    | >90%  |
| 储存条件  | 4-8 °C  |
| 规范性声明 | For research use only. Not for use in diagnostic procedures,                  |

