



CD16 S-pluriBead® anti-hu

Product Images



Description

CD16 pluriBead® suspension for the positive cell isolation of CD16+ target cells from any human sample material. (2ml)

Additional Information

| | |
|-------|---|
| 大小 | 2 ml |
| 品种 | Human |
| 细胞类型 | Neutrophile Granulocytes |
| 样品料 | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| 纯度 | >90% |
| 储存条件 | 4-8 °C |
| 规范性声明 | For research use only. Not for use in diagnostic procedures, |

