




PLT Spin Medium

Product Images



Description

PLT Spin Medium is a ready-to-use, sterile density medium for the isolation of platelets in high yield from fresh human peripheral blood through single step centrifugation.

Additional Information

交付时间	1-2
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
储存条件	Room temperature, protect from light
规范性声明	For research use only. Not for use in diagnostic procedures.
相位范围	6.8 – 7.4

Product Options

Size:	100 ml
	250 ml
	500 ml

