



## ÜberStrainer Set 1

### Product Images



Description

The ÜberStrainer is more than a cell strainer - ÜberStrainer is a sterile Sample Preparation Device to separate, isolate, concentrate and process particle out of complex and/ or diluted suspensions. The air tight screw cap with a Luer Lock connector makes it possible to filtrate an unlimited amount of liquid. The strainer can be easily separated from the rest of the device, fits into a 50 ml tube and into the 20 mm tissue plate cavity.

This set includes 2 pieces of 20 µm, 60 µm, 100 µm, 300 µm, 500 µm.

Additional Information

Mesh Size	Multiple
大小	10 pcs
Fabric material	PET (Polyethylenterephthalat)
外壳材料	PP (Polypropylene)
Mesh Type	woven
不育症	Sterile
运输条件	Room Temperature
储存条件	Room Temperature
规范性声明	For research use only. Not for use in diagnostic procedures.

