



Vortex 1

Product Images



Description

Attractively designed test tube shaker with touch function. Designed for mixing small test samples.

- Small, compact and reliable
- For small containers up to 30 mm in diameter, e.g. test tubes, centrifuge tubes, Eppendorf beakers
- Infinitely adjustable speed from 1000 – 2800 rpm
- Excellent mixing action

- The upper casing and the test tube surface are made from inert plastic. The bottom section is a coated, zinc die-casting

Additional Information

交付时间	1-2
运输条件	Room Temperature
储存条件	Room Temperature

