



## Syringe Strainer 60 µm

### Product Images



Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

Additional Information

交付时间	1-2
Mesh Size	60 µm
大小	12 pcs.
流量控制	Yes
外壳材料	LD-PE (Low Density Polyethylen)
Mesh Type	woven
不育症	Sterile
运输条件	Room Temperature
储存条件	Room Temperature
规范性声明	For research use only. Not for use in diagnostic procedures.

Product Options

Variation:	12 pcs sterile
	12 pcs, non-sterile

