



## M-pluriBead® Maxi Reagent Kit

### Product Images



### Description

Maxi M-pluriBead Reagent Kit contains 100ml 10xWash buffer, stabilization buffer, incubation buffer,

detachment activation buffer, 10x 1,8ml detachment buffer concentrate, 10 M-pluriStrainer, 10 Funnel, 10 connector rings

Additional Information

|                      |                                                                               |
|----------------------|-------------------------------------------------------------------------------|
| Size                 | 10                                                                            |
| Delivery Time (days) | 1-2                                                                           |
| Species              | Human                                                                         |
| Sample Material      | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity               | >90%                                                                          |
| Storage Condition    | Room Temperature; Buffer 4-8 °C; Buffer D -20°C                               |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures,                  |

