



CD326 M-pluriBead® anti-ms

Product Images



Description

CD326 pluriBead® suspension for the positive cell isolation of CD326+ target cells from any murine sample material. (2ml)

Additional Information

Size	2 ml
Species	Mouse (Murine)
Cell Types	Circulating Tumor Cells, CTC, Epithelial cells
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

