



pluriSpin® Human NK Cell Depletion

Product Images



Description

pluriSpin® Human NK Cell Depletion is designed to isolate pure, viable and untouched nk cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels NK cells, while the isolated human nk cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified nk cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

| | |
|----------------------|---|
| Species | Human |
| Cell Types | NK cells |
| Sample Material | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity | >90% |
| Storage Condition | Room Temperature |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

Product Options

| | |
|-------|-------|
| Size: | 2 ml |
| | 10 ml |

