



pluriSpin® Human T Cell Depletion

Product Images



Description

pluriSpin® Human T Cell Depletion is designed to isolate pure, viable and untouched t cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels T cells, while the isolated human t cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified t cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

Delivery Time (days)	1-2
Species	Human
Cell Types	T cells
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

Product Options

Size:	2 ml
	10 ml

