



CD4 M-pluriBead® anti-ms

Product Images



Description

CD4 pluriBead® suspension for the positive cell isolation of CD4+ target cells from any murine sample material. (2ml)

Additional Information

| | |
|----------------------|---|
| Size | 2 ml |
| Species | Mouse (Murine) |
| Cell Types | Monocytes (subset), Helper T cells |
| Sample Material | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity | >90% |
| Storage Condition | 4-8 °C |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

