



CD11b M-pluriBead® anti-ms

Product Images



Description

CD11b pluriBead® suspension for the positive cell isolation of CD11b+ target cells from any murine sample material. (2ml)

Additional Information

Size	2 ml
Species	Mouse (Murine)
Cell Types	Monocytes, Granulocytes, Macrophages, NK cells (subset)
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

