



## CD14 S-pluriBead® anti-ms

### Product Images

---



## Description

---

CD14 pluriBead® suspension for the positive cell isolation of CD14+ target cells from any murine sample material. (2ml)

## Additional Information

---

Size	2 ml
Species	Mouse (Murine)
Cell Types	Monocytes
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

