




PLT Spin Medium

Product Images



Description

PLT Spin Medium is a ready-to-use, sterile density medium for the isolation of platelets in high yield from fresh human peripheral blood through single step centrifugation.

Additional Information

Delivery Time (days)	1-2
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Storage Condition	Room temperature, protect from light
Regulatory Statement	For research use only. Not for use in diagnostic procedures.
pH Range	6.8 – 7.4

Product Options

Size:	100 ml
	250 ml
	500 ml

