




KMO 3 basic (Magnetic Stirrer)

Product Images



Description

The award-winning KMO 3 basic is a small but powerful magnetic stirrer without heating function. It provides:

- unique design with glass surface
- brilliant, clear white display
- illuminated, flat stirring plate
- key-lock button
- reverse rotation
- integrated USB interface (control via labworldsoft and use of Firmware Update Tool)
- recessed thread for support rod
- protection class IP42

Additional Information

Delivery Time (days)	1-2
Shipping Condition	Room Temperature
Storage Condition	Room Temperature

