



ROLLER 6 basic

Product Images



Description

Roller shaker with 6 rolls provides smooth rocking and rolling action at a fixed speed. Used for mixing blood samples, solid and liquid suspensions and viscous samples. The unit can be used with all conventional tubes and cylindrical bottles. Easily removable attachments provide for quick cleaning in case of sample spills.

- Usable at 4 to 50 °C in incubators
- Suitable for continuous operation
- Removable rolls enable shaking of larger tubes
- Easy operation with touch keypad
- Shaking diameter: 24,5 mm (height)
- Speed range: 30 rpm (fixed)

Additional Information

Delivery Time (days)	1-2
Shipping Condition	Room Temperature
Storage Condition	Room Temperature

