



Syringe Strainer 100 μ m

Product Images



Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

Additional Information

Delivery Time (days)	1-2
Mesh Size	100 µm
Size	12 pcs.
Flow Control	Yes
Housing Material	LD-PE (Low Density Polyethylen)
Mesh Type	woven
Sterility	Sterile
Shipping Condition	Room Temperature
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures.

Product Options

Variation:	12 pcs sterile
	12 pcs, non-sterile

