

Mini-Strainer-Aid



Specification

Order No:	43-11000-50
Color:	white
Lot No.:	See label
Material:	Thermoplastic elastomer
Quantity:	50 devices, sterile
Tray material:	Polypropylene
Autoclavable:	yes (121 °C)
Storage:	Room temperature



Component Description

The Mini-Strainer-Aid is designed to facilitate the filtration of small volumes using the Pluriselect Mini-Strainer family. The perforated elastomer aid fits airtight into the Mini-Strainer. The opening holds standard 1 - 5 ml pipetted tips. The elastomer and the specially designed surface create high friction between the pipette, aid, and Mini-Strainer. That allows easy step by step handling. The Strainer-Aid reduces the filtration time in cases where centrifugation is usually needed (small mesh sizes, high suspension density, etc.).

The strainer-aid-pipette-tip-unit can be easily discarded by the tip ejection device of the pipette.

Warning and Limitation

The product is for research and development only, not for diagnostic or therapeutic use.

Pluriselect-usa, Inc. – North America
worldwide
2030 Gillespie Way 103, El Cajon, CA

USA
Phone: +1 619 202 4297
0

support@pluriselect-usa.com

Pluriselect Life Science

Deutscher Platz 5c, D-04103
Leipzig
Germany

Phone: +49 341 333858-

info@pluriselect.com

Pluriselect-usa, Inc. – North America
worldwide
2030 Gillespie Way 103, El Cajon, CA

USA
Phone: +1 619 202 4297
0

support@pluriselect-usa.com

Pluriselect Life Science

Deutscher Platz 5c, D-04103
Leipzig
Germany

Phone: +49 341 333858-

info@pluriselect.com