



Connector Ring

Product Images



Description

The Connector Ring is a sterile device to support active filtration, cell detachment, cell lysis or cell stimulation in combination with a [pluriStrainer](#) and a 50ml standard (conical) centrifuge tube. pluriStrainer, Connector Ring and a 50ml tube can be combined to a unit. This unit allows for a controlled flow through the sieve. When the Luer-lock of the Connector Ring is closed, sample or liquid remains on the sieve. An open Luer-lock starts the separation or the filtration process.

Additional Information

Size	25 pcs.
Delivery Time (days)	1-2
Variation	25 pcs, sterile
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

