

pluriPlix



Order No.:50-01010-80

Description

Universal rotation mixer without own drive system for usage on top of a magnetic stirrer (not included in delivery). Particularly designed for the gentle mixing of cells. Allows for 2 mixing options: over-head and tilting-rotation (plus: adjustable rotation angle).

Additional Information

Size	1 device
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.