



M-pluriBead® Mini Reagent Kit

Product Images



Description

Mini M-pluriBead Reagent Kit contains 100ml 10xWash buffer, stabilization buffer, incubation buffer, detachment activation buffer, 5x 1,8ml detachment buffer concentrate, 5 M-pluriStrainer, 5 Funnel, 5 connector rings

Additional Information

Size	5
Delivery Time (days)	1-2
Species	Human
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	Room Temperature; Buffer 4-8 °C; Buffer D -20°C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

