



universal M-pluriBead® anti-ms

Product Images



Description

pluriBead universal are developed for the use with external antibodies. They are coated with a polyclonal antibody produced from the serum of a host animal. (2ml)

Additional Information

Size	2 ml
Species	Mouse (Murine)
Cell Types	Any
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

