



ÜberStrainer Set 2

Product Images



Description

The ÜberStrainer is more than a cell strainer - ÜberStrainer is a sterile Sample Preparation Device to separate, isolate, concentrate and process particle out of complex and/ or diluted suspensions. The air tight screw cap with a Luer Lock connector makes it possible to filtrate an unlimited amount of liquid. The strainer can be easily separated from the rest of the device, fits into a 50 ml tube and into the 20 mm tissue plate cavity.

This set includes 2 pieces of 1 µm, 5 µm, 10 µm, 30 µm, 60 µm.

Additional Information

Mesh Size	Multiple
Size	10 pcs
Fabric material	PET (Polyethylenterephthalat)
Housing Material	PP (Polypropylene)
Mesh Type	woven
Sterility	Sterile
Shipping Condition	Room Temperature
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures.

