



## Vortex 1

### Product Images



### Description

Attractively designed test tube shaker with touch function. Designed for mixing small test samples.

- Small, compact and reliable
- For small containers up to 30 mm in diameter, e.g. test tubes, centrifuge tubes, Eppendorf beakers
- Infinitely adjustable speed from 1000 – 2800 rpm
- Excellent mixing action

- The upper casing and the test tube surface are made from inert plastic. The bottom section is a coated, zinc die-casting

**Additional Information**

---

Delivery Time (days)	1-2
Shipping Condition	Room Temperature
Storage Condition	Room Temperature

