



## Syringe Strainer 100 $\mu$ m

### Product Images



Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

Additional Information

|                      |                                                              |
|----------------------|--------------------------------------------------------------|
| Delivery Time (days) | 1-2                                                          |
| Mesh Size            | 100 µm                                                       |
| Size                 | 12 pcs.                                                      |
| Flow Control         | Yes                                                          |
| Housing Material     | LD-PE (Low Density Polyethylen)                              |
| Mesh Type            | woven                                                        |
| Sterility            | Sterile                                                      |
| Shipping Condition   | Room Temperature                                             |
| Storage Condition    | Room Temperature                                             |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures. |

Product Options

|            |                     |
|------------|---------------------|
| Variation: | 12 pcs sterile      |
|            | 12 pcs, non-sterile |

