



Multi-Well-Strainer Plates

12 wells, 72 μ m, sterile

Specification

Part No: 43-72123-01
Lot No.: See top bag
Material: Polystyrene
Mesh: Nylon
Mesh surface 1.4 cm² per well
Quantity: 1 device, sterile
Autoclavable: no



Component Description

PluriCell is a product line developed to combine cell propagation with size separation. The multi-well cell nylon-mesh strainers are available in mesh sizes of 72 μ m and 210 μ m. The multi-well strainer plates fit into a wide range of multi-well plates from different manufacturers. They are available for 6- and 12-well plates. Multi-well cell-strainer plates offer a system for culturing and testing of single- or multicell organisms. The strainer plate and the well tray fit together to provide a ready-to-use system for assays measuring paralysis, cytotoxicity, and death of multicellular organisms. Simply transfer parasites, nematodes, or other multicellular organisms into the wells, treat with the target compound, and measure the migration of the organisms.

Warning and Limitation

The product is for research and development only, not for diagnostic or therapeutic use.