




## BSA stock solution (50 ml)

### Product Images



## Description

---

10% BSA stock solution in 1x PBS can used for cell culture, cell separation wash buffer, blocking steps in Western blot, ELISA, IHC and nucleic acid detection methods.

## Features of Blocker BSA

- **Ready to use** - No need to dissolve powder
- **Pure** - High-quality bovine serum albumin as 10% solution provides fewer chances of cross-reaction with assay components than other components
- **Convenient** - Optimal blocking results for specific applications

## Additional Information

---

Size	50 ml
Delivery Time (days)	1-2
Area of interest	Cell Biology
Sterility	Sterile
pH Range	pH 7.4
Shipping Condition	Room Temperature
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

