



pluriSpin® Human CD45⁺ Cell Depletion

Product Images



Description

pluriSpin® Human CD45⁺ Cell Depletion is designed to isolate pure, viable and untouched CD45⁺-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels CD45⁺ cells, while the isolated population free of CD45⁺ are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified CD45⁺ cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

| | |
|----------------------|---|
| Size | 2 ml |
| Species | Human |
| Cell Types | CD45- Cell |
| Storage Condition | 4-8 °C |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures., |

