



Bundle pluriStrainer® Mini 10 µm and Flow Cytometry Tubes

Product Images



Description

Bundle includes sterile pluriStrainers Mini with 10 μm mesh and sterile 5 ml polystyrene (PS) tubes (12x75mm). The combination of a flow cytometry tube and cell strainer offers a convenient way to prepare laboratory samples for flow cytometry analysis. A 10 μm polyethyleneterephthalat (PET) mesh is incorporated into the housing of the pluriStrainer Mini 10 μm . The polystyrene round-bottom tubes are closed with a cap. Tubes and the mini cell strainers will be delivered in separate bags with 25 pcs. each.

Additional Information

Mesh Size	10 µm
Packaging	50 pcs. (25 pcs. / Bag)
Color	turquoise
Flow Control	No
Mesh Fixation	Injected in housing for strong hold
Fabric material	PET (Polyethylenterephthalat)
Housing Material	PP (Polypropylene)
Fits into	1.5 ml reaction tube, 2.0 ml reaction tube, FACS™ tube, Cryo vials 5 ml / 10 ml, reaction tube with lid, 5 ml / 10 ml reaction tube with screw cap, 14- 15 ml conical centrifuge tube, 24 well plate, 48 well plate
Handling	Gripping surface enables easy and aseptic handling.
Sterility	Sterile
Comparable with	Cell Strainer Snap Cap™, Corning Cell Strainer, EASYstrainer™, Cell Sieve, Falcon™ Cell Strainer, Reversible™ Strainer, MACS® SmartStrainer, MACS® Pre-Separation Filters
Stability	RLT® Buffer, TRIzol®, isopropanol, organic solvents
Centrifugable	Yes
Tissue Dissociation	Yes
Disinfectable	Yes
Shipping Condition	Room Temperature
Storage Condition	Room Temperature
Legal information	Buffer RLT® is a trademark of Qiagen, TRIzol® is a trademark of Molecular Research Center Inc., FACS™ is a trademark of BD Biosciences, Corning is a registered trademark of Corning Inc., EASYstrainer™ is a registered trademark of Greiner Bio One International GmbH, Cell Strainer Snap Cap™ and Falcon™ Cell Strainers is a registered trademark of Corning Inc., MACS® SmartStrainer is a registered trademark of Miltenyi Biotec GmbH

