



pluriMag Rack II (15 / 50 ml)

Product Images



Description

pluriMag Rack is a magnetic holder for small-scale protein and nucleic acids purification and cell enrichment or depletion with magnetic beads. pluriMag Rack can be used for sample volumes of up to 50 ml. The holder make an excellent complement to the following pluriMag products (containing magnetic beads) from pluriSelect. pluriMag Rack works also with all standard magnetic micro beads.

It is optimized for efficient magnetic separation of all types of magnetic (micro) beads (range 1–8 μm in diameter) from the supernatant (1 – 45 ml) in seconds, allowing for easy pipette removal without sample loss. The pluriMag Rack combines a strong magnetic field for efficient separations in standard microcentrifuge tubes.

- Best working volume: 1 ml - 45 ml
- Holds 1x 15 ml tube and 1x 50 ml tube
- Good control and visibility of your sample

Additional Information

Size	1 device
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures.,

