



Monocytes Spin Medium, 500 ml

Product Images



Description

Monocytes Spin Medium is a ready-to-use, sterile density medium for the enrichment of monocytes in a high yield from fresh peripheral blood or buffy coat through single step density centrifugation. Depending on the sample a purity of 60-80% is possible. Monocytes Spin is a gentle alternative to bead-mediated monocyte enrichments, especially for downstream applications like cell culture and differentiation.

Additional Information

Size	500 ml
Delivery Time (days)	1-2
Storage Condition	Room temperature, protect from light
Regulatory Statement	For research use only. Not for use in diagnostic procedures.
pH Range	6.8 – 7.4

