




Density Diluent Medium

Product Images



Description

High Density Spin Medium is a ready to use, sterile medium for the preparation of different density working solutions for the enrichment of cells, bacteria and other particles. The working solutions can be prepared by the dilution of High Density Spin Medium with our Density Diluent Medium.

Additional Information

Delivery Time (days)	1-2
Shipping Condition	Room Temperature
Storage Condition	Room temperature, protect from light
Regulatory Statement	For research use only. Not for use in diagnostic procedures.
pH Range	6.8 - 7.4

Product Options

Size:	250 ml
	500 ml

