



pluriStrainer Mini 200 μm (Cell Strainer)

Product Images



Description

The pluriStrainer® Mini is a sieving device (or cell strainer) for a wide range of laboratory applications where filtration and purification of liquids is required. It has been optimized for small sample volumes. Its unique design allows fitting on a wide variation of tubes, e.g. 1.5 ml /2.0 ml reaction tubes, 15 ml conical centrifuge tubes, FACS™ tubes, Cryo vials, 24 well plate and 48 well plate. In combination with a standard flow cytometry tube the pluriStrainer Mini 5µm can be used as replacement for a strainer cap.

Additional Information

Delivery Time (days)	1-2
Mesh Size	200 µm
Packaging	50 pcs. (25 pcs. / Bag) - The tube for flow cytometry is not included.
Color	light transparent brown
Flow Control	No
Mesh Fixation	Injected in housing for strong hold
Fabric material	PET (Polyethylenterephthalat)
Housing Material	PP (Polypropylene)
Mesh Type	woven
Fits into	1.5 ml reaction tube, 2.0 ml reaction tube, FACS™ tube, Cryo vials 5 ml / 10 ml, reaction tube with lid, 5 ml / 10 ml reaction tube with screw cap, 14- 15 ml conical centrifuge tube, 24 well plate, 48 well plate
Handling	Gripping surface enables easy and aseptic handling.
Sterility	Sterile
Comparable with	Cell Strainer Snap Cap™, Corning Cell Strainer, EASYstrainer™, Cell Sieve, Falcon™ Cell Strainer, Reversible™ Strainer, MACS® SmartStrainer, MACS® Pre-Separation Filters
Stability	RLT® Buffer, TRIZOL®, isopropanol, organic solvents
Centrifugable	Yes
Tissue Dissociation	Yes
Disinfectable	Yes
Shipping Condition	Room temperature
Storage Condition	Room temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures

Legal information

Buffer RLT® is a trademark of Qiagen, TRIzol® is a trademark of Molecular Research Center Inc., FACS™ is a trademark of BD Biosciences, Corning is a registered trademark of Corning Inc., EASYstrainer™ is a registered trademark of Greiner Bio One International GmbH, Cell Strainer Snap Cap™ and Falcon™ Cell Strainers is a registered trademark of Corning Inc., MACS® SmartStrainer is a registered trademark of Miltenyi Biotec GmbH

