



Re-Strainer (In-Line-Strainer) 20 μ m - Sterile

Product Images



Description

The Re-Strainer is a filtration device for high volumes of liquid to exclude or concentrate particles from liquids with the possibility of particle recovery. The 6 ml recovery reservoir can be accessed by unscrewing the cap. The Re-Strainer has two female luer-lock-ports.

These ports can be used to connect several Re-Strainer to a filter cascade (male to male luer-lock connector) or integrate the Re-Strainer as an in-line-strainer via luer-lock tube adaptor.

The mesh is an integrated part of the housing to give the device maximum mechanical stability.

Additional Information

Delivery Time (days)	1-2
Mesh Size	20 µm
Packaging	4 devices per bag
Size	20 pcs.
Color	green
Flow Control	yes
Mesh Fixation	Injected in housing for strong hold
Fabric material	PET (Polyethylenterephthalat)
Housing Material	LD-PE (Low Density Polyethylen)
Sterility	Sterile
Stability	Good resistance to non-oxidizing acids and bases, fats, most organic solvents
Shipping Condition	Room Temperature
Storage Condition	Room Temperature
Regulatory Statement	For research use only
Legal information	The product is for research and development only, not for diagnostic or therapeutic use

