



## SnapCap Steel Extraction Strainer, non-sterile, 10 pcs

### Product Images

---



### Description

---

The SnapCap Extraction Strainer is screwed into a standard 50 ml polypropylene tube. All liquids for which PP and Nylon are resistant can be used as extraction medium.

## Additional Information

---

Delivery Time (days)	1-2
Mesh Size	1 $\mu\text{m}$
Size	10 pcs
Fabric material	Non-magnetic SS 316 steel
Housing Material	LD-PE (Low Density Polyethylen)
Sterility	Sterile
Storage Condition	Room Temperature
Lot No. / EXP	on bag, top
Metal Ring	Magnetic steel SS304
Strainer volume	9 cm <sup>3</sup>
Cap sealing	Thermoplastic elastomer 60A
50 ml tube	PP
Snap cap material	Nylon 6
Glue stability	Sensitive to Toluene and Acetone
Legal information	The product is for research and development only, not for diagnostic or therapeutic use.

## Product Options

---

Variation:	50 $\mu\text{m}$
	500 $\mu\text{m}$

