



## PluriStrainer-Maxi (Bottle Top), 5 pcs. , single packed, non-sterile

### Product Images

---



### Description

---

The pluriStrainer maxi can connect to GL32 or GL80 bottles via our adaptor. The various mesh sizes and the

unique stackability of the pluriStrainer maxi allows to build a filter cascade and separate sample material into different fractionations. Furthermore, the pluriStrainer maxi has a special lid for an automated refilling. The sample material can be drawn directly from the sample vessel into the strainer via this lid. A tubing port enables low pressure for filtration and automated refilling.

### Additional Information

---

Delivery Time (days)	1-2
Size	5
Housing Material	LD-PE (Low Density Polyethylen)
Sterility	Non-Sterile
Storage Condition	Room Temperature

### Product Options

---

Variation:	5 µm
	20 µm
	40 µm
	70 µm
	100 µm
	200 µm
	300 µm
	400 µm
	500 µm
	750 µm
	1000 µm
	1400 µm
	2000 µm

